



# TRIPLE THERAPY

Le sorgenti in dotazione al sistema **TRIPLE THERAPY** consentono la più ampia scelta ai medici e terapisti, i quali potranno utilizzare un'unica apparecchiatura come **YAG-laser**, utilizzando esclusivamente la sorgente con lunghezza d'onda a 1064 nm o come **laser di elevata potenza** utilizzando la sorgente da 10W con lunghezza d'onda 810 nm. Il sistema **TRIPLE THERAPY** ottimizza gli spazi all'interno dello studio medico e ottimizza l'investimento economico in un'unica apparecchiatura (anziché due apparecchi con 2 lunghezze d'onda differenti).

Gli studi effettuati sono stati condotti utilizzando contemporaneamente sia lo YAG laser a 1064 nm che il diodo a 810 nm con differenti frequenze. I risultati sono stati rapidi ed efficaci.

La rotazione della testata consente agli operatori di trattare le parti più difficili del corpo umano (spalla, zona cervicale della colonna vertebrale) senza costringere il paziente a posizioni scomode.

La possibilità di effettuare una termografia sulla zona da trattare consente di valutare i valori termici nella fase antecedente la somministrazione della TRIPLE THERAPY e verificare subito dopo le variazioni termiche indotte dalla terapia.



*The laser sources supplied with **TRIPLE THERAPY** give a wide possibility of choice to Doctors and Therapists, which can be free to use in a single device a **YAG laser** (working with wavelength 1064 nm) and a **HIGH POWER laser** (working with wavelength 810 nm – 10W). Providing two different wavelengths in a single device, the **TRIPLE THERAPY** system is a real medical office space and economical investment optimizer.*

*The studies carried out using both 1064 nm YAG laser source and 810 nm diode with different frequencies have shown quick and effective results.*

*The rotating head facilitate the operators while treating the most complicated parts of the body (shoulder, cervical area along the vertical spine) without uncomfortable postures for the patient.*

*The possibility to make a pre-treatment thermography on the area to heal and after the TRIPLE THERAPY application, allow to verify the temperature variations, resulted by therapy.*

**level**  
Medical Systems

# TRIPLE THERAPY



**level**

Level s.r.l.

Via Mattarella, 59 - 30037 Gardigiano di Scorzè (VE)  
T. +39 041 449132 - level@level-laser.it - www.level-laser.com

UNI EN ISO 9001-2008  
UNI CEI EN ISO 13485-2012  
DIRECTIVE 93/42/EEC - Annex II

extremonline.it

DP25 05/15

TRIPLE THERAPY

**level**

Camera termografica con software dedicato  
Thermographic camera with specific software

